

Product Datasheet

Usutu Virus (USUV-NS1) Antibody HM623

Product Name	Usutu Virus (USUV-NS1) Antibody
Catalog Number	HM623
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Antigen (Vienna 2001)
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No cross-reactivity with other Flaviviruses
Specificity	Usutu Virus
Matched Pairs	HM624 (detection)
Preservatives	0.05% Sodium Azide