

Product Datasheet

Usutu Virus (USUV-NS1) Antibody HM624

| | |
|----------------------------|--|
| Product Name | Usutu Virus (USUV-NS1) Antibody |
| Catalog Number | HM624 |
| Description | Ion Exchange Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, LF |
| Immunogen | Recombinant Antigen (Vienna 2001) |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | > 90 % |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Cross Reactivity | No cross-reactivity with other Flaviviruses |
| Specificity | Usutu Virus |
| Matched Pairs | HM622 or HM623 (capture) |
| Preservatives | 0.05% Sodium Azide |