

Product Datasheet

Fentanyl Antibody HM628

| | |
|----------------------------|---|
| Product Name | Fentanyl Antibody |
| Catalog Number | HM628 |
| Description | Ion Exchange Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, LF |
| Immunogen | Fentanyl-KLH |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | > 90 % |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Cross Reactivity | Acetyl Fentanyl 100%; Norfentanyl & Acetyl Norfentanyl 90%; Alfentanil 30%; Risperidone 20%; Carfentanil 5% |
| Specificity | Fentanyl |
| Matched Pairs | LA395 (capture) |
| Preservatives | 0.05% Sodium Azide |