

Product Datasheet

Helicobacter pylori Antibody HM629

Product Name	Helicobacter pylori Antibody
Catalog Number	HM629
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF, TIA
Immunogen	Recombinant H. pylori (LA396)
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Campylobacter, Clostridium difficile, E. coli, L. monocytogenes, Salmonella, Shigella, Y. enterocolitica
Specificity	H. pylori
Matched Pairs	HM629 (self-pairs)
Preservatives	0.09% Sodium Azide