

Product Datasheet

Rotavirus (VP6) Antibody HM630

Product Name	Rotavirus (VP6) Antibody
Catalog Number	HM630
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Rotavirus VP6 (LA397)
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short term: +2-8 °C Long term: -20°C
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Campylobacter, Clostridium difficile, E. coli, L. monocytogenes, Salmonella, Shigella, Y. enterocolitica, Enterovirus, Norovirus, Astrovirus, Adenovirus, Cryosporidium parvum, Giardia lamblia, H. pylori, Entamoeba hystolitica, S. aureus
Specificity	Rotavirus VP6
Matched Pairs	HM630 (self-pairs)
Preservatives	0.09% Sodium Azide