

Product Datasheet

Clostridium difficile Glutamate Dehydrogenase (GDH) Antibody HM631

Product Name	Clostridium difficile Glutamate Dehydrogenase (GDH) Antibody
Catalog Number	HM631
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant C. difficile GDH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Campylobacter, E. coli, H. pylori, L. monocytogenes, Salmonella, Shigella, Y. enterocolitica, S. aureus
Specificity	C. difficile GDH
Matched Pairs	HM631 (self-pairs)
Preservatives	0.09% Sodium Azide