

## Product Datasheet

### Clostridium difficile Glutamate Dehydrogenase (GDH) Antibody HM631

<b>Product Name</b>	Clostridium difficile Glutamate Dehydrogenase (GDH) Antibody
<b>Catalog Number</b>	HM631
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF
<b>Immunogen</b>	Recombinant C. difficile GDH
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	+2-8 °C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95 % by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No reactivity with Campylobacter, E. coli, H. pylori, L. monocytogenes, Salmonella, Shigella, Y. enterocolitica, S. aureus
<b>Specificity</b>	C. difficile GDH
<b>Matched Pairs</b>	HM631 (self-pairs)
<b>Preservatives</b>	0.09% Sodium Azide