

Product Datasheet

Norovirus genogroup II Antibody HM638

Product Name	Norovirus genogroup II Antibody
Catalog Number	HM638
Species/Host	Mouse
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Group II capsid protein
Matched Pairs	HM885 (detection)
Preservatives	0.05% Sodium Azide