

Product Datasheet

Human MICA ELISA Kit 850.910.096

Product Name	Human MICA ELISA Kit
Catalog Number	850.910.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	156.25-5000 pg/ml
Sensitivity	123 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with IL-7, IL-12p40, IL-2R, Trail, CD31, IL-13Rα2, IL-5, gp130, TNFα.
Specificity	Recognizes both natural and recombinant human MICA
Sample Types	Supernatant, Buffered Solutions. Not suitable for use with serum and plasma.