

Product Datasheet

Entamoeba histolytica Antibody HM643

Product Name	Entamoeba histolytica Antibody
Catalog Number	HM643
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	E. histolytica
Matched Pairs	HM643 (self-pairs)
Preservatives	0.1% Sodium Azide