

Product Datasheet

Estrone-3-Glucuronide (E1G) Antibody HM658

| | |
|----------------------------|---|
| Product Name | Estrone-3-Glucuronide (E1G) Antibody |
| Catalog Number | HM658 |
| Species/Host | Rat |
| Tested Applications | LF, ELISA |
| Immunogen | Estrone-3-Glucuronide-KLH |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | >90% by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl |
| Cross Reactivity | Estriol <6%; Estriol-3-glucuronide, Estradiol-3-glucuronide <5%; 17 α -Ethinylestradiol <3%; 17 α -hydroxyprogesterone, Testosterone, Progesterone, Pregnenediol-3-glucuronide, Aldosterone, Androstenedione, Estrone sulfate <1% |
| Specificity | Estrone and Estrone-3-glucuronide |
| Matched Pairs | LA432 (capture) |
| Preservatives | 0.05% Sodium Azide |