

Product Datasheet

Human IL-23 ELISA Kit 850.920.096

Product Name	Human IL-23 ELISA Kit
Catalog Number	850.920.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	156.25-5000 pg/ml
Sensitivity	20 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant human IL-23
Sample Types	Serum, Plasma, Cell culture supernatant