

Product Datasheet

Enterovirus 71 Antibody HM673

Product Name	Enterovirus 71 Antibody
Catalog Number	HM673
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with Enterovirus 70, Poliovirus 2, Coxsackie A, Coxsackie B
Specificity	Enterovirus 71 VP1
Preservatives	0.1% Sodium Azide