

Product Datasheet

C-Reactive Protein, Canine (cCRP) Antibody HM686

Product Name	C-Reactive Protein, Canine (cCRP) Antibody
Catalog Number	HM686
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Canine CRP
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Canine CRP
Matched Pairs	HM685 (detection)
Preservatives	0.02% Proclin 300