

# Product Datasheet

## Myoglobin (Mb) Antibody HM688

|                     |                                                                           |
|---------------------|---------------------------------------------------------------------------|
| Product Name        | Myoglobin (Mb) Antibody                                                   |
| Catalog Number      | HM688                                                                     |
| Species/Host        | Mouse                                                                     |
| Tested Applications | LF, ELISA, CLIA                                                           |
| Immunogen           | Recombinant Human Myoglobin                                               |
| Brand               | EastCoast Bio                                                             |
| Form/Appearance     | Protein A Purified Monoclonal Antibody                                    |
| Storage             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype             | IgG                                                                       |
| Clonality           | Monoclonal                                                                |
| Purity              | >95% by SDS-PAGE                                                          |
| Buffer              | 10 mM Phosphate Buffered Saline, pH 7.4                                   |
| Specificity         | Human Myoglobin                                                           |
| Matched Pairs       | HM687 (capture)                                                           |
| Preservatives       | 0.02% Proclin 300                                                         |