

Product Datasheet

D-dimer Antibody HM694

Product Name	D-dimer Antibody
Catalog Number	HM694
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, CLIA, LF
Immunogen	Recombinant D-dimer
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Human D-dimer
Matched Pairs	HM695 (detection)
Preservatives	0.02% Proclin 300