

# Product Datasheet

## Human CD62P ELISA Kit 850.930.096

Product Name	Human CD62P ELISA Kit
Catalog Number	850.930.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	4.4-140 ng/ml
Sensitivity	3.6 ng/ml
Tested Applications	ELISA
Alternative Names	P-selectin
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant human P-Selectin
Sample Types	Serum, Plasma, Cell culture supernatant, Amniotic fluid