

Product Datasheet

ST2 Antibody HM696

| | |
|----------------------------|--|
| Product Name | ST2 Antibody |
| Catalog Number | HM696 |
| Description | Protein A Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA, CLIA |
| Immunogen | Recombinant ST2 |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | > 95 % by SDS-PAGE |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Specificity | Human ST2 |
| Matched Pairs | HM697 (detection) |
| Preservatives | 0.05% Proclin 300 |