

Product Datasheet

Procalcitonin (PCT) Antibody HM699

Product Name	Procalcitonin (PCT) Antibody
Catalog Number	HM699
Species/Host	Mouse
Tested Applications	LF, ELISA, CLIA
Immunogen	Recombinant procalcitonin
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Human procalcitonin
Matched Pairs	HM698 (capture)
Preservatives	0.02% Proclin 300