

Product Datasheet

Myoglobin (Mb) Antibody HM705

Product Name	Myoglobin (Mb) Antibody
Catalog Number	HM705
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified human myoglobin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Human Myoglobin
Matched Pairs	HM704 (capture)
Preservatives	0.09% Sodium Azide