

Product Datasheet

Norovirus genogroup II Antibody HM712

Product Name	Norovirus genogroup II Antibody
Catalog Number	HM712
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Norovirus GII.4 VLP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Norovirus GII
Matched Pairs	HM713 (detection)
Preservatives	0.09% Sodium Azide