

## Product Datasheet

### Norovirus genogroup II Antibody HM715

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Norovirus genogroup II Antibody  |
| <b>Catalog Number</b>      | HM715  |
| <b>Description</b>         | Purified Monoclonal Antibody   |
| <b>Species/Host</b>        | Mouse  |
| <b>Tested Applications</b> | ELISA, LF  |
| <b>Immunogen</b>           | Norovirus GII.4 VLP  |
| <b>Brand</b>               | EastCoast Bio  |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody   |
| <b>Storage</b>             | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| <b>Isotype</b>             | IgG1   |
| <b>Clonality</b>           | Monoclonal   |
| <b>Purity</b>              | > 95 % by SDS-PAGE   |
| <b>Buffer</b>              | Phosphate Buffered Saline, pH 7.0  |
| <b>Specificity</b>         | Norovirus GII  |
| <b>Matched Pairs</b>       | HM714 (capture)  |
| <b>Preservatives</b>       | 0.09% Sodium Azide   |