

Product Datasheet

Serum Amyloid A (SAA) Antibody HM718

Product Name	Serum Amyloid A (SAA) Antibody
Catalog Number	HM718
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant Serum Amyloid A
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Human Serum Amyloid A
Matched Pairs	HM719 or HM720 (detection)
Preservatives	0.09% Sodium Azide