

Product Datasheet

Human IL-17A ELISA Kit 850.940.192

Product Name	Human IL-17A ELISA Kit
Catalog Number	850.940.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	3.125-1000 pg/ml
Sensitivity	3 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reaction with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-17A
Sample Types	Serum, Plasma, Cell culture supernatant