

Product Datasheet

Heparin Binding Protein (HBP) Antibody HM721

Product Name	Heparin Binding Protein (HBP) Antibody
Catalog Number	HM721
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA, CLIA
Immunogen	Purified Human Heparin Binding Protein
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Heparin Binding Protein
Matched Pairs	HM722 (capture or detection)
Preservatives	0.05% Sodium Azide