

Product Datasheet

Heparin Binding Protein (HBP) Antibody HM721

| | |
|----------------------------|---|
| Product Name | Heparin Binding Protein (HBP) Antibody |
| Catalog Number | HM721 |
| Description | Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, CLIA, LF |
| Immunogen | Purified Human Heparin Binding Protein |
| Alternative Names | Azurocidin |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | +2-8 °C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | > 95 % by SDS-PAGE |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Specificity | Heparin Binding Protein |
| Matched Pairs | HM722 (capture or detection) |
| Preservatives | 0.05% Sodium Azide |