

Product Datasheet

Heparin Binding Protein (HBP) Antibody HM722

Product Name	Heparin Binding Protein (HBP) Antibody
Catalog Number	HM722
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, CLIA, LF
Immunogen	Purified Human Heparin Binding Protein
Alternative Names	Azurocidin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short-Term: 2-8 °C; Long-Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Heparin Binding Protein
Matched Pairs	HM721 (capture or detection)
Preservatives	0.05% Sodium Azide