

Product Datasheet

Streptococcus Group B Antibody HM731

Product Name	Streptococcus Group B Antibody
Catalog Number	HM731
Species/Host	Rabbit
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No cross-reactivity with other Streptococcal groups
Specificity	Group B Streptococcus species
Matched Pairs	HM731 (self-pairs)
Preservatives	0.1% Sodium Azide