

Product Datasheet

Enterobacteriaceae Antibody HM732

Product Name	Enterobacteriaceae Antibody
Catalog Number	HM732
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	OMP from Shigella species mixture
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with Pseudomonas aeruginosa, Vibrio cholera and L. monocytogenes
Specificity	Enterobacterial Common Antigen (ECA)
Preservatives	0.1% Sodium Azide