

Product Datasheet

Anti-Mullerian Hormone (AMH) Antibody HM733

Product Name	Anti-Mullerian Hormone (AMH) Antibody
Catalog Number	HM733
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, CLIA
Immunogen	AMH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	Hemoglobin, LH <2%; Inhibin B <0.5%
Specificity	AMH
Matched Pairs	HM737 (detection)
Preservatives	0.1% Proclin 300