

Product Datasheet

Anti-Mullerian Hormone (AMH) Antibody HM740

Product Name	Anti-Mullerian Hormone (AMH) Antibody
Catalog Number	HM740
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, CLIA, LF
Immunogen	Recombinant AMH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	> 95 %
Buffer	Phosphate Buffered Saline, pH 7.0
Specificity	AMH
Matched Pairs	HM1050 (detection)
Preservatives	0.09% Sodium Azide