

## Product Datasheet

### C-Reactive Protein (CRP) Antibody HM751

<b>Product Name</b>	C-Reactive Protein (CRP) Antibody
<b>Catalog Number</b>	HM751
<b>Description</b>	Protein A Purified Monoclonal Antibody
<b>Immunogen</b>	Purified human CRP
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified monoclonal antibody
<b>Storage</b>	Short term: 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95% by SDS-PAGE
<b>Buffer</b>	Phosphate-buffered saline, pH 7.0 ± 0.2, 0.09% sodium azide as a preservative.
<b>Specificity</b>	Human CRP
<b>Matched Pairs</b>	HM752 (detection)