

Product Datasheet

Hemoglobin (Hb) Antibody HM753

| | |
|----------------------------|--|
| Product Name | Hemoglobin (Hb) Antibody |
| Catalog Number | HM753 |
| Description | Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Human hemoglobin |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | > 95 % by SDS-PAGE |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.0 |
| Specificity | Human hemoglobin |
| Matched Pairs | HM061 or HM782 (detection) |
| Preservatives | 0.09% Sodium Azide |