

Product Datasheet

Hemoglobin (Hb) Antibody HM753

Product Name	Hemoglobin (Hb) Antibody
Catalog Number	HM753
Description	Purified Monoclonal Antibody
Immunogen	Purified human hemoglobin
Brand	EastCoast Bio
Form/Appearance	Purified monoclonal antibody
Storage	Short term: 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	≥ 95 % by SDS-PAGE
Buffer	Phosphate-buffered saline, pH 7.0 ± 0.2, 0.09% sodium azide as a preservative.
Specificity	Human hemoglobin
Matched Pairs	HM061 or HM782 (detection)