

Product Datasheet

Epstein-Barr Virus (EBV-EBNA-1)-HRP Antibody HM757

Product Name	Epstein-Barr Virus (EBV-EBNA-1)-HRP Antibody
Catalog Number	HM757
Species/Host	Goat
Tested Applications	ELISA, WB, IHC
Immunogen	Purified recombinant EBNA-1
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody-HRP Conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 with 10 mg/ml BSA
Specificity	EBNA-1
Preservatives	0.002% Thimerosal