

# Product Datasheet

## Human IP-10 ELISA Kit 850.950.192

Product Name	Human IP-10 ELISA Kit
Catalog Number	850.950.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	6.25-200 pg/ml
Sensitivity	5.7 pg/ml
Tested Applications	ELISA
Alternative Names	CXCL10
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reaction with other human cytokines
Specificity	Recognizes both natural and recombinant human IP-10
Sample Types	Serum, Plasma, Cell culture supernatant