

Product Datasheet

Nipah Virus (NiV-gF) Antibody HM768

Product Name	Nipah Virus (NiV-gF) Antibody
Catalog Number	HM768
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Nipah Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.15M Potassium Phosphate, 0.85% NaCl, pH 7.2
Matched Pairs	HM766 or HM769 (capture)
Preservatives	0.05% Sodium Azide