

Product Datasheet

Toxoplasma gondii (p30) Antibody HM772

Product Name	Toxoplasma gondii (p30) Antibody
Catalog Number	HM772
Species/Host	Mouse
Tested Applications	ELISA, WB, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	T. gondii p30 (surface antigen 1)
Preservatives	0.1% Sodium Azide