

Product Datasheet

Gastrin 17 (G-17) Antibody HM775

Product Name	Gastrin 17 (G-17) Antibody
Catalog Number	HM775
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA, CLIA, IFA
Immunogen	Recombinant Gastrin 17
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	LOD 0.3 pmol/L in IFA
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Gastrin 17
Matched Pairs	HM774 (capture)
Preservatives	0.02% Proclin 300