

# Product Datasheet

## Human IL-27 ELISA Kit 850.960.192

|                     |  |
|---------------------|--|
| Product Name        | Human IL-27 ELISA Kit  |
| Catalog Number      | 850.960.192  |
| Description         | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range               | 31.25-1000 pg/ml   |
| Sensitivity         | 12.8 pg/ml   |
| Tested Applications | ELISA  |
| Alternative Names   | NULL   |
| Brand               | Diaclone   |
| Storage             | Store at +4°C  |
| Cross Reactivity    | No cross reactivity with other human cytokines   |
| Specificity         | Recognizes both natural and recombinant human IL-27  |
| Sample Types        | Serum, Plasma, Cell culture supernatant  |