

Product Datasheet

Feline Immunodeficiency Virus (FIV-p24) Antibody HM781

Product Name	Feline Immunodeficiency Virus (FIV-p24) Antibody
Catalog Number	HM781
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, IHC, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Note	Addition of 0.2% Triton-X recommended in non-reducing applications such as ELISA
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with FeLV, Feline Herpesvirus type 1, Feline Coronavirus, Feline Calicivirus
Specificity	FIV p24 (mature protein and higher MW precursors)
Matched Pairs	HM780 (capture)
Preservatives	0.05% Sodium Azide