

Product Datasheet

Progastrin-releasing Peptide (ProGRP) Antibody HM784

Product Name	Progastrin-releasing Peptide (ProGRP) Antibody
Catalog Number	HM784
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, CLIA, CMIA
Immunogen	Recombinant ProGRP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffer with 150 mM NaCl, pH 7.4
Matched Pairs	HM785 or HM788 (detection)
Preservatives	0.1% Proclin 300