

Product Datasheet

Progastrin-releasing Peptide (ProGRP)-HRP Antibody HM788

Product Name	Progastrin-releasing Peptide (ProGRP)-HRP Antibody
Catalog Number	HM788
Description	Purified Monoclonal Antibody-HRP Conjugate
Species/Host	Mouse
Tested Applications	ELISA, CLIA, CMIA
Immunogen	Recombinant ProGRP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody-HRP Conjugate
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffer, pH 7.4 with 10% NBCS and 50% glycerol
Matched Pairs	HM784 (capture)
Preservatives	0.1% Proclin 300