

Product Datasheet

Mayaro Virus (MAYV-E1) Antibody HM798

| | |
|----------------------------|---|
| Product Name | Mayaro Virus (MAYV-E1) Antibody |
| Catalog Number | HM798 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Recombinant Mayaro Virus |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | >90% by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride |
| Matched Pairs | HM795, HM797 or HM798 (capture) |
| Preservatives | 0.05% Sodium Azide |