

Product Datasheet

Mayaro Virus (MAYV-E1) Antibody HM798

Product Name	Mayaro Virus (MAYV-E1) Antibody
Catalog Number	HM798
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Mayaro Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM795, HM797 or HM798 (capture)
Preservatives	0.05% Sodium Azide