

Product Datasheet

Astrovirus Antibody HM803

Product Name	Astrovirus Antibody
Catalog Number	HM803
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Astrovirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No cross reactivity with Adenovirus, C. difficile, Cryptosporidium, Enterovirus, E. coli, Giardia, H. pylori, L. monocytogenes, Norovirus, Rotavirus, Y. enterocolitica
Matched Pairs	HM804 (capture or detection)
Preservatives	0.09% Sodium Azide