

Product Datasheet

Campylobacter Antibody HM805

Product Name	Campylobacter Antibody
Catalog Number	HM805
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG2b
Clonality	Monoclonal
Purity	90 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with <i>H. pylori</i> , <i>E. coli</i> , <i>S. enteritidis</i> , <i>S. paratyphi</i> , <i>S. typhi</i> , <i>S. typhimurium</i> , <i>S. boydii</i> , <i>S. dysenteriae</i> , <i>S. flexneri</i> , <i>S. sonnei</i> , <i>Y. enterocolitica</i> , <i>L. monocytogenes</i> , <i>S. aureus</i>
Matched Pairs	HM806 (capture)
Preservatives	0.09% Sodium Azide