

## Product Datasheet

### Campylobacter Antibody HM806

<b>Product Name</b>	Campylobacter Antibody
<b>Catalog Number</b>	HM806
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95 % by SDS-PAGE
<b>Buffer</b>	Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No reactivity with <i>H. pylori</i> , <i>E. coli</i> , <i>S. enteritidis</i> , <i>S. paratyphi</i> , <i>S. typhi</i> , <i>S. typhimurium</i> , <i>S. boydii</i> , <i>S. dysenteriae</i> , <i>S. flexneri</i> , <i>S. sonnei</i> , <i>Y. enterocolitica</i> , <i>L. monocytogenes</i> , <i>S. aureus</i>
<b>Matched Pairs</b>	HM805 (detection)
<b>Preservatives</b>	0.09% Sodium Azide