

Product Datasheet

Cryptosporidium Antibody HM807

Product Name	Cryptosporidium Antibody
Catalog Number	HM807
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified Native <i>C. parvum</i>
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with <i>E. dispar</i> , <i>E. histolytica</i> , <i>G. lamblia</i>
Specificity	Cryptosporidium
Matched Pairs	HM807 (self-pairs)
Preservatives	0.09% Sodium Azide