

Product Datasheet

Lactoferrin Antibody HM810

Product Name	Lactoferrin Antibody
Catalog Number	HM810
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with bovine lactoferrin, human/pig/bovine hemoglobin, human calprotectin, human/pig/bovine transferrin
Matched Pairs	HM811 (capture)
Preservatives	0.09% Sodium Azide